

Leaders In GC Sample Introduction Technology







Purge and Trap Sample Introduction Instruments

Reliable Purge and Trap Concentrator with 2nd Generation Autosampler Options for both Water and Soil Samples

465 Limestone Road, Oxford, PA 19363-0277, USA Tel: 800-541-6593 Fax: 610-932-4158 www.cdsanalytical.com



Introduction

If gas chromatography plays a role in your daily operations, you already know CDS Analytical as the world's leading manufacturer of pyrolysis instrumentation. But we also manufacture worldclass Purge & Trap GC front-end sampling instruments.

The CDS-7000E was acknowledged as a workhorse in the industry for decades as the most cost effective Purge and Trap concentrator. With the introducion of the 2nd generation autosampler 7350 and 7450 families, this dependable Purge and Trap concentrator is unlocking its full productivity while maintaining the unrivalled chromatographic resolution and the industrial leading Cost Performance Index.



Figure 1: VOC's at 20 ppb concentration purged from water. Figure 2: Multiple runs of the six gases at 500 ppt. Figure 3: Linearity for Vinyl Chloride from 500 ppt to 50 ppb

	Disblaradiflueromathana		
	Dichlorodifluoromethane Chloromethane	32.	1,
	Vinyl chloride	33.	CI
	Bromomethane	34.	
		35.	
	Chloroethane	36.	
	Trichlorofluoromethane	37.	
	1,1-Dichloroethylene	38.	0-
	Methylene chloride	39.	St
	trans-1,2-Dichloroethene		Tr
10.	1,1-Dichloroethane	41.	Is
	2,2-Dichloropropane	42.	Br
12.	cis-1,2-Dichloroethene	43.	
13.	Bromochloromethane	44.	
14.	Chloroform		
	1,1,1-Trichloroethane	46.	
16.	Carbon tetrachloride	47.	
17.	1,2-Dichloropropene		
	Benzene		
	1,2-Dichloroethane		
20.	Flurorbenzene (I.S.)		
21.	Trichloroethylene	52.	
22.	1,2-Dichloropropane		
23.	Dibromomethane	54.	
24.	Bromodichloromethane		
	cis-1,3-Dichloropropene		
26.	Toluene	57.	
27.	trans-1,3-Dichloropropene		
28.	1,1,2-Trichloroethane		
29.	Tetrachloroethylene		He
30.			
	Dibromochloromethane		

- lorobenzene ,1,2-Tetrachloroethane
- hylbenzene -Xylene Xylene Xylene

- 1,2,2-Tetrachloroethane 2,3-Trichloropropane

- chlorotoluene 3,5-Trimethyl benzene Chlorotoluene rt-Butyl benzene

- rt-Butyl benzene 24-Trimethyl benzene 3-Dichlorobenzene Isopropyl toluene 4-Dichlorobenzene 4-Dichlorobenzene, d-4 2-Dichlorobenzene 2-Dichlorobenzene
- - loroporpane 3.5-Trichlorobenzene

 - ohthalene .3-Trichlorobenzene

7000E Design Highlights

- The sample path is made of SilcoNert 2000 coated stainless steel, which is the best chemically inert material commercially available
- Valve oven reaches as high as 350°C, which is capable of covering the widest range of analytes including Gases, VOCs, low boilers, medium boilers and semi volatiles.
- Hydroguard coated Wet Trap follows EPA method and eliminates over 96% of moisture
- Patented foam sensor prevents possible contaminations
- Quick access panel makes regular maintenance easier than ever
- Front panel shows the instrument status at a glance
- Powerfull Windows based user interface makes it compatible with leading GC systems



7000E User Interface



7000E Specifications

- Programmable times: 0– 999.9 minutes
- Trap: 0.3 cm OD x 28.5 cm Length
- Transfer Line: 1.5 m Long, 0.020" ID Silconert tubing in flexible heated jacket
- Desorb pre-heat option
- Maximum Operating Temperatures
 - Valve oven 350°C
 - Transfer Line 350°C
 - Wet Trap 300°C
- Adsorbent Trap 300°C
- Vocarb or Tenax Adsorbent Trap
- Weight: 8.6 kg
- 25 cm W x 44 cm H x 46 cm D

Options:

- Hot water rinse module
- Electronic mass flow controller

2nd Generation CDS Autosampler Platform

The new generation CDS autosampler platform promises a more precise handling for sample loading. This 2nd generation platform will unlock a new era of GC front-end sampling instruments with the following features:



(1) Improved Precision: The position accuracy of the 2nd generation autosampler is improved by 50%. Due to a new high precision drive shaft and components, 3D (X, Y and Z) placement accuracy is now within 1 mm as compared to 1.5 mm in the first generation. (2) Quicker Operation: The new autosampler platform does not need to HOME after each movement of the vial or tube. There is less movement needed with the pick-and-place arm which ultimately should extend it's useful life.

(3) Ease of Alignment: The alignment calibration has shifted from hardware calibration to software calibration. This greatly improves and accelerates installation set up.

(4) Longer Life Expectancy. The new control board has adopted a SMART acceleration curve to lower the impact on the driving mechanism by 40%. This extends the life of the system as less strain will be applied on the motors and drive shaft.



osampler												
mmunication												
omm: COM1 ~	Rack1: 24 V Custom			Water sample position			Waste position					
communication					х	15		х	0			
	Rack2:		~	Custom	Y	0	Test	×	0		Test	
	Rack3:	24	~	Custom		·	Update		-	- 10	Iodate	
option	Rack4: 24		Ý	Custom	Z	50	0,0300	z	50		opone	
innotation:	Rack			first postion			current position					
1. click up before move to another position	۲	1 02	03	04	x	180			x	Y Y	z	
to another position					Y	518	Test					
2.Reference position:					*	310	Undate					
Rack1:					z	165	Update					
First Pos:180-518-165	Accor	ding to	the Coo	dinate mobile						single step		
End Pos:924-2978-165	х	0							+	5		
Rack2:	¥	0		Test		lest a	ostion					
First Pos: 1294-508-165		Ľ.,			×	919	coben		10			
					x	919	Test			Up Dow	n	
End Pos:2039-2973-165		Up		Down	Y	2978						
Reck3:						Update						
First Pos:2424-508-165					z	Z 165						
End Pos:3169-2958-165			Home									
010 105:3103-2330-103						save						

Sample Vial Position Calibration through the Software

CDS-7350 Purge and Trap Water Autosampler









Make CDS Your Universal Inlet Partner



 CDS Analytical Headquarters:

 465 Limestone Road P.O. Box 277 Oxford, PA 19363-0277
 Tel: 800.541.6593
 610.932.3636
 www.cdsanalytical.com